

## KS-TY9000 Optical Emission Spectrometers



(Picture just for reference)

### 1. Summary

KS-TY9000 is Full Range of Solutions for the Entire Metals Industry. It use full-digital technology to replace bulky photo multiplier tube (PMT) simulation technology and keep pace with international spectrometer technology. The adoption of vacuum optical chamber design, full-digital excitation light source, advanced CCD detectors, and high-speed data readout system equips the device with high properties, ultra-low limit of detection (LOD), long-term stability and repeatability. The analytical precision can meet the requirements of laboratory standard, the analytical date is stable and reliable. Widely used in Metallurgy, casting, machinery processing and other industries incoming and outgoing product quality control.

### 2. Scope of supply

#### 2.1 Supply list

No.	Item	Specifications	QTY	Unit	Remarks
1	Optical Emission Spectrometer	KS-TY9000	1	Set	Include in quotation.
2	Business computer		1	Set	
3	Printer		1	Set	

4	Tungsten electrode		1	Pc
5	Electrode pressure spring		1	Pc
6	The electrode fixed screw		2	Pcs
7	Electrode brush	Φ6	2	Pcs
8	Lens holder sealing ring	Φ16×2.65	2	Pcs
9	Lens holder sealing ring	Φ72×2.65	1	Pc
10	Exhaust gas filter core		1	Pc
11	Degreasing cotton		15	G
12	Power socket		1	Pc
13	M2 Internal hexagonal wrench		1	Set
14	M4 Internal hexagonal wrench		1	Set
15	Long handle phillips screwdriver		1	Pc
16	13 – 15 fixed wrench		1	Pc
17	Spark machine screw	M6	4	Pcs
18	Fuse	10A	2	Pcs
19	Argon gas pressure reducing valve		1	Pc
20	Argon gas pipeline		2	M
21	Exhaust pipe		1	Pc
22	Gas bottle		2	Pcs
23	Print paper		1	Pc
24	Spectrometer manual		1	Set

Notes: The Buyer knows and agrees that the Seller may adjust the supply when the production of the above-mentioned instruments, spare parts and other suppliers change due to technological progress, product upgrading, market changes, etc.

**2.2 Optional items:**

Optional item	Qty	Note
High purity argon gas (pure $\geq$ 99.999%)	1 bottle	Should prepare by customer. If can't get should use argon purifier to get 99.999% argon gas.
Small air conditioning	1 set	Essential, customer should prepare by themselves.
High precision magnetic saturation voltage stabilizer (1KVA)	1 set	Need for voltage instability, if the customer uses intermediate frequency furnace, then need 3KVA Voltage Stabilizer
Argon purifier	1 set	When argon gas is not as pure as 99.999%.
Samples grinder (ferrous metals) or Mini lathe (non-ferrous metals)	1 set	Essential
Type samples	pieces	Essential for foundry

Please select according to your requests:

### OES Analytical Program

No	Elements	Low Alloy Steel	Cast Iron
1	C	0.006-1.6	1.81-4.5
2	Si	0.0014-2.57	0.15-3.96
3	Mn	0.021-19.74	0.114-3.42
4	P	0.0032-0.135	0.021-0.82
5	S	0.0013-0.111	0.003-0.211
6	Cr	0.002-12.42	0.03-2.5
7	Ni	0.0107-4.0	0.01-4.77
8	Mo	0.004-3.98	0.003-1.7
9	Al	0.003-1.29	0.002-0.214
10	Cu	0.0032-0.7	0.06-1.73
11	Co	0.0088-0.5	0.0046-0.14
12	Ti	0.001-0.297	0.004-0.5
13	Nb	0.0057-0.506	0.003-0.33
14	V	0.0082-1.26	0.0014-1.13
15	W	0.05-1.97	0.001-0.97
16	Mg		0.0001-0.1
17	B	0.0017-0.0083	0.002-0.26
18	Sn	0.0035-0.054	0.0007-0.563
19	Zn	0.002-0.04	0.001-0.0055
20	As	0.0019-0.056	0.002-0.26
21	Bi	0.0002-0.013	0.0017-0.05
22	Zr	0.005-0.41	
23	La		0.0002-0.023
24	Ca	0.0001-0.055	
25	Ce	0.0002-0.019	0.0002-1.81
26	Sb	0.0004-0.01	0.004-0.2
27	Fe	Reference	Reference

**Notes:**

1. The elements highlighted in orange are available normally, those in blue can be added according to the standard sample available or not.
2. The analysis curves highlighted in green can be added according to current stand samples in the markets, those highlighted in orange are available normally, and it can meet 80% customer's need.
3. Numbers of Elements added or channels are depended on the price.
4. For special material beyond the detection range of elements and measurement, buyer and supplier can make a decision on the internal control standard.
5. Due to the lag of the national standard material, the elements will be carried out according the original standard before new material coming out.







6. The final explanation right of element table and measuring range belongs to KASON Group.






### 3. Technical Data of KS-TY9000

#### 3.1 Parameter

	Item	Index	
Optical System	Focal Length	400mm	
	Wavelength range	165~589nm (Extendable)	
	Detector	High resolution CCD Multi detectors	
	Degree of vacuum	Auto control within 6-20 pa	
	Pixel resolution	30pm	
	Grating line	2400m1/mm	
	First order spectral line dispersion rare	1.2nm/mm	
	Average resolution ratio	10pm/pixel	
	Full spectrum		
	Light room temperature is controlled automatically		
Spark Source	Type	Digital arc and spark source	
	Spark frequency	100-1000HZ	
	Discharge current	1-400A	
	Ignition voltage	>15000V	
	Excitation light	Optimization of discharge parameters design	
		High energy precombustion technology HEPS	
Processor	High-speed data synchronization acquisition and processing		
Spark Stand	Electrode	Tungsten electrode technology	
	Make up	Thermal deformation self-compensation design	
	Argon flushed with minimal consumption of Argon		
	Spray discharge electrode technology		
	Adjustable electrode technology		
Others	Measurable elements	Fe、Al、Cu、Ni、Ti、Co、Zn、Sn、Mg、Pb etc	
	Dimension	800mm(L)*700mm*470mm(H)	
	weight	About 100kg	
	Storage temperature	0℃-45℃	
	Operating temperature	10℃-35℃, 23±2℃ is recommended	
	Power	AC220V/50Hz(Customized)	
	Power consumption	Excitation:700W/Stand by:100W	
	Argon quality	99.999%, Argon pressure>4Mpa	
	Argon consumption	5L/min during spark mode	
	Interface	Ethernet data transmission based on DM9000A	

#### 3.2 Main features

<p><b>1. Optimized vacuum optics system</b></p> <p>(1) Integrated optics room and Paschen Runge construction design, making all the spectrum lines focused on the gratings.</p> <p>(2) Direct-jet type optics technology and MgF2 material lens to make sure ultraviolet wave's best energy of elements, such as C, S, P and N.</p>	
<p><b>2. Automatic Light-Path Correction</b></p> <p>(1) With automatic light path correction, optical system automatically scans the spectral lines to ensure the correctness of received lines and avoid tedious scanning of wave peaks.</p> <p>(2) The instrument automatically identifies specific spectral lines and compares them with original stored lines to determine the location of the drift and find the present pixel position for analysis among the lines.</p>	
<p><b>3. Single-Board Lens Design</b></p> <p>(1) The adopted specialized entrance window separating from vacuum in the vacuum optical system can be operated under the system working status. The adopted single-board lens structure in the optical lens is convenient for routine cleaning and maintenance.</p> <p>(2) It is not required to maintain the device in the daily operation, and there is no consumable and renewal part.</p>	
<p><b>4. Optical Chamber Integration</b></p> <p>(1) Specialized optical chamber structural design makes the volume of the chamber smaller, with less than half air exhaust speed of ordinary spectrometers.</p> <p>(2) Integrated design and high-precision processing in the vacuum chamber improves the duration of vacuum.</p>	
<p><b>5. Vacuum Anti-Oil-Returning Technology</b></p> <p>(1) Multi-level separating vacuum anti-oil-returning technology absorbs vacuum compaction and baffle valves to ensure the complete separation of vacuum optical chamber from vacuum pump during non-operation time.</p> <p>(2) The intermediate addition of vacuum oil filtering device ensures the oil inside the vacuum pump not to enter the vacuum chamber, and ensures CCD detectors and optical components to work under reliable situation.</p>	
<p><b>6. Open-Access Inspire Stand</b></p> <p>(1) Flexible sample clamp design of Open-Access Excitation Stand satisfies user on-site analysis of samples with different sizes and shapes</p> <p>(2) The small sample clamps in cooperative use can make the analytic precision of wire rods achieve 1.5 mm at the minimum.</p>	

<p><b>7. Injection Electrode Technology</b></p> <p>(1) The instrument adopts the most internationally-advanced injection electrode technology with tungsten used. Under the status of excitation, the electrode forms argon gas injection flow surrounded. Thus, there is no opportunity for the surrounding excitation points to contact external air so as to improve precision of excitation.</p> <p>(2) The attached specialized argon gas channel design significantly reduces the usage amount of argon gas and use cost for customers.</p>	
<p><b>8. Integrated Gas Channel Block</b></p> <p>(1) The excitation stand is made of alloys with good heat dissipation to achieve solidity, durability and cleaning convenience.</p> <p>(2) The gas supply system adopts integrated gas channel block and electrode self-flushing function to create good environment for excitation.</p>	
<p><b>9. Full Digital Inspire Light Source</b></p> <p>(1) The system uses the most internationally-advanced plasma inspire light source, and generates ultra-stable energy release to excite samples in the environment filled with argon gas.</p> <p>(2) Full Digital Inspire Light Source ensures ultra-high resolution and high-stability output rate of plasma in the excitation samples.</p> <p>(3) Full Digital Inspire Light Source can satisfy inspire requirements for different type of materials.</p>	
<p><b>10. High-Speed Data Acquisition</b></p> <p>(1) The instrument adopts high-performance CCD devices, UV coating technology and high-performance FPGA, DSP and ARM processors.</p> <p>(2) The system has ultra-high functions of data acquisition and analysis, and can realize automatic real-time monitoring and control of the block operation status of optical chamber temperature, vacuum degree, argon gas pressure, light source and excitation chamber.</p>	
<p><b>11. Ethernet Data Transfer</b></p> <p>(1) Ethernet cards and TCP / IP protocol are connected between computers and spectrometers to avoid electromagnetic interference, fiber aging defects. Meanwhile, the computers and printers are completely placed in external position for the benefit of promotion and substitution.</p> <p>(2) Complete network systems.</p> <p>(3) The system can remotely monitor the status of devices, operate the system through multiple channels, as well as</p>	



## 19 YEARS

Professional focused on testing equipment


KASON is established in 2003, owns more than 8000 square meters factory, has a professional sales teams, modern enterprise technology center, scientific and technological research and development team.

Machines passed the European CE authentication, American FDA certificate and ISO 9001.

Products sold to USA, Canada, Australia, Europe, Africa etc, more than 130 countries and supply OEM service for many customers

## PROFESSIONAL TEAM

KASON has a professional sales teams, modern enterprise technology center, scientific and technological research and development team.



## OUR CUSTOMERS (more than 130 countries)



