

Electric Marking Machine for Tensile Test



(Picture just for reference)

Product Description

The motorised rebar marking machine is designed for steel reinforcement elongation marking instruments. It is more accurate than manual marking. The efficiency can be increased by several times. This machine has features of precise height, easy operation and convenient movement. It is essential instrument for the laboratory.

Standard: ISO 6935-1, ISO 6935-2

Main Features:

* Easy clamping for sample mounting

- * Marking length of 300 mm in increment of 10 mm
- * Up to 36mm diameter rebar

Specification

Product Name	Electric steel bar marking machine for rebar tensile test
Model	KASON-DQ5-10
Marking Length	300mm
Print Gauge Length	5mm/10 mm
Specimen Diameter	6-36mm
Dimension	830x 250x300mm
Weight	45KG
Power	220V



Authentication Certificate of Quality Management System

Registration No.19922Q01072R0S

This is to certify that

Jinan Kason Testing Equipment Co., Ltd.

Unified Social Credit Code: 9137010458994939X0

Registered/Production Address: No. 4715, Jingshi West Road,
Huaiyin District, Jinan City, Shandong Province

Office Address: Room 2-1108, Harmony Square, Huaiyin District, Jinan City

Is in conformity with :

GB/T 19001-2016/ISO 9001:2015

The certificate is valid for the following scope :

Assembly and sales of testing machines

This certificate should be used effectively with the legal requirements of the administrative license, qualification license, compulsory certification, which the certification scope involved.

In the case that the organization regularly receives surveillance assessments, the certificate shall be valid when used together with the Quick Response Code.

Date of Issue: 2022-07-19

Date of Expiry: 2025-07-18

Date of Initial Certification: 2022-07-19

General Manager



中国认可
国际互认
管理体系
MANAGEMENT SYSTEM
CNAS C182-M

China Standard Huaxin (Beijing)
Certification Center Co.,Ltd.



Address of Certification Authority: 5 Floor, YuQuan Building, No.3 Shijingshan Road, Shijingshan District, Beijing Tel: 010-88255986 P.C.: 100049

Note: The detailed information is available on the Certification and Accreditation Administration of the People's Republic of China official website (www.cnca.gov.cn) the authority's website(www.cshcc.cn) and the Quick Response Code on the left.



MACHINERY DIRECTIVE ATTESTATION OF CONFORMITY

Technical file of the company mentioned below has been inspected and audit has been completed successfully.

2006/42/EC Machinery Directive Annex VIII has been taken as references for these processes.

Company Name	: Jinan Kason Testing Equipment Co., Ltd.
Company Address	: No.4715, Jingshi West Road, Huaiyin District, Jinan, Shandong, China
Related Directives and Annex	: 2006/42/EC Machinery Directive/Annex VIII
Related Standards	: EN ISO 12100:2010; EN 60204-1:2018
Product Name	: Testing Machine
Report No and Date	: TCF-JKTE230605-MD
Product Brand/Model/Type	: WDS-1, WDS-2, WDS-3, WDS-5, WDS-10, WDS-20, WDS-30, WDS-50, WDS-100, WDS-200, WDS-300, WDW-1E, WDW-2E, WDW-3E, WDW-5E, WDW-10E, WDW-20E, WDW-30E, WDW-50E, WDW-100E, WDW-200E, WDW-300E, WDW-500E, WDW-600E, WDW-1000E, ETM103, ETM203, ETM303, ETM503, ETM104, ETM204, ETM304, ETM504, ETM105, WE-300D, WE-600D, WE-1000D, WE-2000D, HUT305S, HUT605S, HUT105S, HUT126S, WEW-300D, WEW-600D, WEW-1000D, WEW-1200D, WEW-2000D, WEW-3000D, HUT305E, HUT605E, HUT105E, HUT126E, HUT-A306E, HUT-A206E, WAW-300D, WAW-600D, WAW-1000D, WAW-1200E, WAW-2000E, WAW-3000, WAW-2000A, WAW-3000A, WAW-4000A, HUT-A206W, HUT-A306W, HUT-A406W, HUT305W, HUT605W, HUT105W, HUT126W, WAL-LAW, HUT-L, JB-300, JB-300B, JB-500B, JBS-300, JBS-500, JBW-300, JBW-500, JBDW-300, JBDW-500, CST-50, VU-2Y, VU-1S, DWC-40, DWC-60, DWC-80, DWC-100, DWC-196, YES-1000, YES-2000, YES-3000, YAW-1000, YAW-300, YAW-2000, YAW-3000, YAW-5000, NDS-200, NDS-500, NDS-1000, NDS-2000, NDW-200, NDW-500, NDW-1000, NDW-2000, NDW-3000, NDW-5000, NDW-10000, EZ-1, EZ-3, EZ-10, JWJ-10, TLS, TLW, HWS, PLS-20, PLS-50, PLS-100, PLS-300, PLS-500, PLS-1000, KASON-1, RDJ-10, RDJ-20, RDJ-50, RDJ-100, DWTT-20000, DWTT-30000, DWTT-40000, DWTT-50000, DWTT-80000, DWTT-100000, PTF-302, PTF-302M, PTF-502, PTF-502M, PTF-302T, PTF-502T, PTF-302W, PTF-502W, PTF-502TW, PTF302CW, PTF452CW, PTF602CW, PTF752CW, PTF302CT, PTF452CT, PTF602CT, PTF752CT, PTF302W, PTF452W, PTF602W, PTF752W, PTF302CW, PTF452CW, PTF602CW, PTF752CW, PTF302CT, PTF452CT, PTF602CT, PTF752CT, HCT106, HCT206, HCT306, KSCIM, HYJ, HEZ-3, HEZ-5, HEZ-10, HEZ-20, KS8T-10, KASON 370.02, KASON 370.05, KASON 370.10, KASON 370.25, KASON 370.50, MMRW, GBW-60B, GBS-60B, HUT, ETM, XNR, GW-40, GW-50, MRS, MMU, MRH, GBW, GW, LSW, HBT, GLY, GCJ-1000, GWQ-1000, GNZ-1000, GQR-1000, GJR-300, GSS-1000, GMS-100, GTW-500, GQC-5, GGG-100, HR-150A, HRD-150, HRS-150, XHR, HB, HVS, HV, HBRV, W-20, HBA-1, LHT, Q-ZA, MP-2, XQ-2B, XD-30M, 4XCE, 4XB, JSZ, MOPAO, MPT, P, MP, BG, LBG-32, MPJ-35, Q-Z, Q-100B, SQ, Q-3A, QG-4, LSQ, LDQ, GTQ-5000A, GTQ-5000B, DTQ, ZXQ, 4XC, HST, KASON, KS, KSCUT, HTR, HTPR, HTB, HSTR-45S, HTV, HTMV, HTBRV, METRPRESS, MPJ

Certificate Number	: M.2023.206.C86313
Initial Assessment Date	: 07.06.2023
Registration Date	: 08.06.2023
Reissue Date/No	: -
Expiry Date	: 07.06.2028

UDEM International Certification
Auditing Training Centre Industry
and Trade Inc. Co.

The validity of the certificate can be checked through www.udem.com.tr. Upon completion of EC declaration of conformity it is used solely at the manufacturer's responsibility. This certificate remains the property of UDEM International Certification Auditing Training Centre Industry and Trade Inc. Co. to whom it must be returned upon request. The above named firm must keep a copy of this certificate for 15 years from the registration of certificate. This certificate only covers the product(s) stated above and UDEM must be noticed in case of any changes on the product(s)

Address: Muthukent Mahallesi 2073 Sokak (Eski 93 Sokak) No:10 Çankaya - Ankara - TURKEY

Phone: +90 0312 443 03 90 **Fax:** +90 0312 443 03 76

E-mail: info@udem.com.tr www.udem.com.tr



FOCUS IN MATERIAL TEST

KASONTTEST®

JINAN KASON TESTING
EQUIPMENT Co., LTD.

Duandian Industrial Park , Jingshi Road, Jinan City, China.

P: +86 159 10081986

E: admin@jnkason.com | W: www.syjlab.com

