

# KASON M-1

## Metallographic Specimen Pre-grinding Machine

GRINDING AND POLISHING MACHINE



### ★ TECHNICAL PARAMETERS

Model	KASON M-1
Pre-grinding disc	Ø230mm
Number of grinding and polishing discs	Single disk
Rotating speed	900rpm
Motor	90W
Dimensions and Net weight	380*490*350mm,11kg
Package size and Gross weight	575*465*540mm,25kg
Voltage	220V, 50Hz, Single phase

### ★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Metallographic sample polishing machine	1	Water inlet pipe	1
Dust cover	1	Product Manual	1
Metallographic sand paper Ø230mm	2	Product certification	1
Drain pipe	1		



Product dimensions

Product details



Metallographic consumables



### ★ PRODUCT FEATURES

- Speed control system optimized for grinding and polishing machines
- Easy-to-operate switches and knobs for easy operation
- Finely ground and surface-treated work disc to ensure smooth sample surface
- The integrally formed high-strength fiberglass material shell is sturdy and durable and will never rust.