

# KASON-MPT

## Semi-automatic Metallographic Grinding Polishing Head

GRINDING AND POLISHING MACHINE



### ★ TECHNICAL PARAMETERS

Model	KASON-MPT
Rotate Speed	50 rpm
Sample Force	0~40N
Sample Capacity	1~3 PIECES
Dimension and Weight	183*209*252mm, 13kg
Packing size and Gross weight	400*300*550mm, 18kg
Power Supply	Single phase, 220V, 50Hz or three phases, 380V, 50HZ



Metallographic products

Product details



Metallographic consumables



### ★ PRODUCT FEATURES

- MPT polishing and grinding head is designed according to investigation and researches on the market and customers' requirements. MPT can prepare one, two or three samples at one time. It can be fixed onto many models of polishing and grinding machines produced by us, It is easy-to-use, and the quality of finished sample is high. MPT is an ideal option for the lab of factories, science and research institutions and universities, which can prepare sample by unmanned operation and increase the efficiency.