

KASON-AKT Abrasion Tester For Rubber & Shoe Sole Rubber Wear Resistance tester



Product instruction

Applied to test the wear ability of vulcanized rubber. It tests the abrasion consumption in certain testing range by the friction under certain tilt angle and load. With digital setting and displaying the times of friction, automatically stop and integrate the host and controller.

Standards

ISO:4661-1:1993, ISO 2481:1988, ISO 471:1983, ISO 1826:1981

Features:

1. Digital setting and display of abrasion times, auto stop
2. Nice appearance, easy to operate

Specifications:

Model	KASON-AKT324A	KASON-AKT324B
Grinding Wheel	Outside dia Φ 150mm,25mm(w),grit 36#	Outside dia Φ 150mm,25mm(w),grit 40#
Grinding wheel	34 \pm 1r/min	35 \pm 1r/min
Specimen Rotating	76 \pm 2r/min	1kg. 6lb
Specimen Slanted	0~35° adjustable	
Counter	LCD, 0~99,999,999	
Power	220V \pm 10% 50HZ, Or specified	

FOCUS IN MATERIAL TEST

KASONTTEST®

**JINAN KASON TESTING
EQUIPMENT Co, LTD.**

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | **W:** www.syjlab.com

