

KASON-5020P Pneumatic Rubber Die Cutter Tensile Sample Making Machine



(Picture just for reference)

19 YEARS

Development Base of Advanced Testing Machine Of China

Product description:

KASON-5020P pneumatic punching press machine using pneumatic power to cut a variety of test strips, suitable for laboratory cutting rubber, plastic, leather and so on. Test dumbbell test pieces, and other special materials, can also act as a small press machine.

Specifications:

Model	KASON-5020P
Types	pneumatic
Moving stroke	20-60mm
Air pressure	0.6MPa
Working platform	170mmX150mm
Dimension	310*470*690mm
Net weight	60kg

KASON reserves the right to modify the technical and stetics characteristics included in this document, without previous notice

DuandianIndustrial Park ,
Jingshi Road, Jinan City,
Shandong Province, China.

Tel. +86 531 58595086
Fax. +86 531 86113769
Email: info@jnkason.com