

KASON-TG10 Rubber Thickness Gauge



(Picture just for reference)

Application:

Rubber Thickness Gauge is mainly used for measuring the thickness of vulcanized rubber and plastic products. It has the following characteristics: easy structure, easy usage, reliable quality and removable base.

Standards:

BS-903-A38, ISO-4648, AASTMD1813, ISO 2589, SATRA TM48

Specification:

Measure the range: 0-10mm, resolution value is 0.01mm;

Have the diameter of examining the head: \varnothing 2mm; \varnothing 5mm; \varnothing 10mm; \varnothing 4mm*1mm

Dimension: L*W*H(mm)138*60*180;

Weight quality: 130g, 110g, 10g

Weight: about 1.5kg

19 YEARS

Professional focused on testing equipment


KASON is established in 2003, owns more than 8000 square meters factory. has a professional sales teams, modern enter prise technology center, scientific and technological research and development team.

Machines passed the European CE authentication, American FDA certificate and and ISO 9001.

Products sold to USA, Canada, Australia, Europe, Africa etc, more than 130 countries and supply OEM service for many customers

PROFESSIONAL TEAM

KASON has a professional sales teams, modern enter prise technology center, scientific and technological research and development team.



OUR CUSTOMERS (more than 130 countries)



