

# Soft Candy Hardness Tester



**19 YEARS**

Development Base of Advanced Testing Machine Of China

## Product introduction

It can be used as a food brittleness tester and can be used to test the elasticity of samples such as jelly, fish balls, fish cakes, fish mince, minced meat and meat balls.

try. Can be used to test the hardness and brittleness of small snacks such as twists, potato chips, and gummies.

## Feature

1. Using chip control and high-precision load cell, high accuracy.
2. Use high-performance motor and ball screw drive to ensure the transmission of the instrument.
3. Large touch screen with English display and friendly human-machine interface operation. Complete testing fully automatically. Easy to operate, with test data statistics office Management function, micro printer output.
4. Special probes and accessories.

## Technical Parameters

Measuring range: (0~250)N or (0~25)kgf

Resolution: 0.01N

Indication accuracy:  $\leq \pm 1\%$ ,

Displacement accuracy: 0.01mm;

Shearing speed: 0.1~20 mm/s;

Adjustable cylindrical probe:  $12.7 \pm 0.1$  mm;

Sample bottle specifications:  $59 \pm 1$  mm (inner diameter) \* 85 mm (height);

Voltage: 220V 50Hz

## 19 YEARS

Professional focused on testing equipment

KASON is established in 2003,owns more than 8000 square meters factory.has a professional sales teams, modern enter prise technology center,scientific and technological research and development team.

Machines passed the European CE authentication,American FDA certificate and and ISO 9001.

Products sold to USA, Canada, Australia, Europe, Africa etc,more than 130 countries and supply OEM service for many customers


## PROFESSIONAL TEAM

KASON has a professional sales teams, modern enter prise technology center,scientific and technological research and development team.





## METALLGURGICAL PRODUCT SHOW



*KASON reserves the right to modify the technical and stetics characteristics included in this document, without previous notice*

DuandianIndustrial Park ,  
Jingshi Road, Jinan City,  
Shandong Province, China.

Tel. +86 531 58595086  
Fax. +86 531 86113769  
Email: admin@jnkason.com