

**KASON HTW-1500
Wood Hardness Tester**



5mm
carbide ball indenter

PRODUCT FEATURES:

- Set the pressing depth through the panel: 5.64 or 2.82 mm
- Set the loading speed through the panel: 3, 4, 5, 6 mm/min
- Automatically display the hardness value

TECHNICAL PARAMETERS:

Model	KASON HTW-1500
Maximum load	15,000N
Load reading accuracy	1N(One order of magnitude higher than the standard requirement (standard requirement: 10N))
Indentation depth accuracy	1microns(One order of magnitude higher than the standard requirement (standard requirement:0.01mm))
Power supply	AC220V/300W
Size	550*210*750mm
weight (Kg)	89

STANDARD CONFIGURATION:

Item	Qty	Item	Qty
indenter	1 piece	Instrument box	1 piece
Dustproof plastic cover	1 piece	power cable	1 piece
Instruction Manual	1 piece	Product certification	1 piece