

KASON M-300 Abrasive Belt Grinding Machine



Product instruction

- 1.Easy to operate, use weak current, safe and convenient.
- 2.Fully closed shell, easy to clean, motor isolation, prevent water and debris erosion.
- 3.With cooling device, can prevent dust damage to human body and workpiece burns.
- 4.Stable operation, automatic correction of sand belt, not easy to run off-track.

Main Specifications

Item	Specifications
Center distance of drum	310 mm
Abrasive belt	Perimeter: 920 mm, width: 100 mm, standard with 60# grit
Surface velocity of sand belt	790 m/min
Motor	Power: 550 W, Speed: 1400 r/min
Power supply	380 V, 50 Hz
Cooling device	Tap water cooling
Overall dimensions	690 × 520 × 355 mm
Weight	49 kg

FOCUS IN MATERIAL TEST

KASONTTEST®

**JINAN KASON TESTING
EQUIPMENT Co, LTD.**

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | **W:** www.syjlab.com

