

# METPRESS-2A

## Automatic Metallographic Specimen Mounting Press

METALLOGRAPHIC MOUNTING PRESS



### ★ PRODUCT FEATURES

- Electric loading of pressure
- Equipped with water cooling system
- Touch screen LCD operation

### ★ TECHNICAL PARAMETERS

Model	METPRESS-2A
Mould diameter	22mm or 30mm or 45mm(Choose One )
Holding time range	0-99min
Heating temperature range	35-200°C
Cooling method	Water cooling
Total input power	950W
Dimension and net weight	485*457*435mm,51.8kg
Package Size and Gross Weight	610*546*687mm,64.5kg
Power supply	Single phase, 220V, 50Hz

### ★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Metallographic Specimen Mounting Press	1Set	Pipe hoop	2pc
Plastic spoon	1PC	Manual	1pc
Thermocouple	1PC	Certification	1pc
Water inlet pipe, Water outlet pipe	1PC		



### Metallographic products



### Metallographic sample hot Inlaid powder



### Metallographic consumables