

METPRESS-3A

Automatic Metallographic Specimen Mounting Press



★ TECHNICAL PARAMETERS

Model	METPRESS-3A
The Sample Diameter	Ø22mm, Ø30mm, Ø40mm, Ø45mm (choose one, Two samples at a time)
Heater	1500W
Total power consumption	1700W
Heating temperature range	35~200
Cooling method	Water Cooling
Dimension and Net Weight	395*460*560mm,50kg
Package Size and Gross Weight	610*495*670mm,64kg
Power supply	single phase, 220V, 50Hz

★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Metallographic Specimen Mounting Press	1Set	Thermocouple	1pc
Plastic spoon	1PC	Water inlet pipe, Water outlet pipe	Each 1pc
Thermocouple	1PC	Manual	1pc
Pipe hoop	2PC	Certification	1pc



Metallographic products

Product details



Metallographic consumables



★ PRODUCT FEATURES

- Specially designed mould locking cap
- Touch screen operation, automatic mounting process
- Equipped with water cooling system