

METPRESS-5A

Automatic Metallographic Specimen Mounting Press

METALLOGRAPHIC MOUNTING PRESS



★ TECHNICAL PARAMETERS

Model	METPRESS-5A
The Sample Diameter	Ø22mm or Ø25mm or Ø30mm or Ø40mm or Ø45mm(choose 1)
Sample preparation pressure	0-10MPa
Sample preparing time	5-6min
Temperature adjusting range	10-300 °C
Cooling method	Water Cooling
Data storage	10 groups
Heater	1600W (22mm, 30mm), 2000W (45mm)
Dimension and Net Weight	435*610*430mm,52kgs
Package Size and Gross Weight	715*530*585mm,64kgs
Power Supply	220V, 50Hz,Single phase

★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Metallographic Specimen Mounting Press	1Set	Pipe hoop	2pc
Plastic spoon	1PC	Filter	1pc
Thermocouple	1PC	Certification	1pc
Water inlet pipe, Water outlet pipe	Each 1pc	Manual	1pc



Product details

Metallographic products



Metallographic consumables



★ PRODUCT FEATURES

- Quick sample preparation within 5 minutes
- 10 groups of data storage
- Electro-hydraulic pressure loading method