

METPRESS-6P

Full Automatic Metallographic Specimen Pneumatic Mounting Press

METALLOGRAPHIC MOUNTING PRESS



★ PRODUCT FEATURES

- Simple structure, quick sample preparation
- Pneumatic pressure applying & releasing
- Automatic sample preparing, automatic sample pushing out

★ TECHNICAL PARAMETERS

Model	METPRESS-6P
Mould diameter	Ø22mm, or Ø30mm, or Ø45mm (choose one mould)
Heater	1200W
Sample preparing time	10-15 minutes
Temperature setting range	0-200°C
Cooling temperature range	0-999min
Constant temperature time	0-99min
Air pressure requirement	0.4-0.7Mpa
Total power supply	1000W
Cooling method	Water cooling
Dimension and Net Weight	435×520×605mm,77kg
Power supply	single phase, 220V, 50Hz

★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Metallographic Specimen Mounting Press	1Set	Pipe hoop	2pc
Plastic spoon	1PC	Filter	1pc
Thermocouple	1PC	Certification	1pc
Water inlet pipe, Water outlet pipe	Each 1pc	Manual	1pc



Metallographic products



Metallographic consumables



Metallographic consumables