

# METPRESS-PVM1

## Vacuum Mounting Press

METALLOGRAPHIC MOUNTING PRESS



### ★ TECHNICAL PARAMETERS

Model	KASON METPRESS-PVM1
Vacuum chamber size	Φ260*125 mm
Vacuum	up to 85 kPa
Sample mold diameter	25 ~ 50 mm
Travel speed	5 rpm
Power	6 W
Dimensions	403*460*400 mm
Power supply	220V, 50Hz, Single phase

### ★ PRODUCT FEATURES

- METPRESS-PVM1 is an innovative solution that is used in high-tech laboratories and production centers. They allow you to easily and safely carry out the process of a sample under pressure. This model has a wide range of capabilities due to its functionality. The diameter of the vacuum chamber is d. 260 mm, which guarantees its large capacity for various types of samples. The speed of movement is 5 rpm, which means that you will get a quick and effective result in the operation of the device. There are no difficulties in operating this product due to its reliable design and intuitive interface.



### Metallographic products



### Metallographic cold inlay



### Metallographic consumables