

KASON MPJ1B Specimen Abrasive Wheel Grinding Machine**Product instruction**

KASON-MPJ1B spectrum grinding machine is mainly used for surface preparation of spectroscopic analysis samples and metal, ceramics, glass, etc. The whole machine has a compact structure, single disc and single speed, and can be used for grinding with grinding wheels and sandpaper of different particle sizes and hardness; polishing is the preparation A main process in the metallographic sample process, the machine rotates smoothly, with low noise, convenient operation and high work efficiency. It is suitable for metallographic laboratories of factories and scientific research units.

Main Specifications

| | |
|------------------------|---------------|
| Model | KASON-MPJ1B |
| Grinding disc diameter | Φ400mm |
| Speed | 1400r/min |
| Input voltage | 2.2KW |
| Motor | AC380V 50HZ |
| Dimension | 650×650×950mm |
| Net weight | 93kg |

FOCUS IN MATERIAL TEST

KASONTTEST®

JINAN KASON TESTING
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | **W:** www.syjlab.com

