

KASON-PS-20 Fully Automatic Gypsum Hardness Tester



19 YEARS

Development Base of Advanced Testing Machine Of China

Product introduction

The hardness of gypsum was measured according to the method specified in GB/T17669.3-1999 Determination of mechanical properties of building gypsum and ISO 3051 standard.

This hardness tester has high precision, fully automatic measurement, and can get the hardness value with just one click.

Feature

- High measurement accuracy: depth measurement reaches 0.2 microns, load accuracy: 0.2%.
- Controlled by computer; fully automatic, hardness results can be obtained with one click.
- Automatically eliminate unqualified indentation depth values.
- Customized user reports.

Technical parameters

Depth accuracy	0.2 μ m
Load accuracy	0.20%
Hardness resolution	0.01 N/mm ²
Measuring range	1-6000
size	520*240*720(mm)
power supply	AC220V/100W
weight	45kg

Standard configuration

host	1 set	Test bench	1 piece
indenter	1 set	10mm hard steel ball	10 pieces
brand computer	1 set	Hardness measurement software	1 set
RS-232 communication line	1 piece	power cable	1 piece
Instrument box	1 set	Dustproof plastic cover	1 piece
Instruction Manual	1 piece	Product certification	1 piece

19 YEARS

Professional focused on testing equipment

KASON is established in 2003,owns more than 8000 square meters factory.has a professional sales teams, modern enter prise technology center,scientific and technological research and development team.

Machines passed the European CE authentication,American FDA certificate and and ISO 9001.

Products sold to USA, Canada, Australia, Europe, Africa etc,more than 130 countries and supply OEM service for many customers


PROFESSIONAL TEAM

KASON has a professional sales teams, modern enter prise technology center,scientific and technological research and development team.





METALLURGICAL PRODUCT SHOW



KASON reserves the right to modify the technical and stetics characteristics included in this document, without previous notice

DuandianIndustrial Park ,
Jingshi Road, Jinan City,
Shandong Province, China.

Tel. +86 531 58595086
Fax. +86 531 86113769
Email: admin@jnkason.com