

Moisture Content Testing Machine



Product instruction

Moisture Content Testing Machine using the principle of heating and drying, moisture values measured object is automatically displayed, the data is stable and accurate.

Measurement: powder form, granular, paste, liquid, organic and inorganic chemical raw materials, a variety of cereals, grain, seed, feed, starch, food, meat, tobacco, tea, paper, grass, soil, sludge, textiles, fuel, resin, silica, plastic particles, coal, fertilizer, medicine, ceramics, wood, ore.

Technical data:

Model	KASON-110MW	KASON-105MW	KASON-102MW	KASON-100MW	KASON-200MW
Range(g)	110	110	110	110	200g
Readability(g)	0.01	0.005	0.002	0.001	0.001g
Operating	5°C---35°C				
Display	LCD				
Calibration	External				
Pan size (mm)	Φ90				
Heat source	Halogen lamp				
Temp.sensor	PT-100				
Temp.range	40°C---199°C				
Temp.step	1°C				
Moisture range	0.0%-100%		0.00%---100%		
Dry residual range	100.00%---0.00%				
ATRO	100%---999%				
Time setting	1---99min by 1s step				
Autostop	0.1%---9.9%within 1 min				
Nuber of storage	15				
Interface	R232				

FOCUS IN MATERIAL TEST

KASONTTEST®

JINAN KASON TESTING
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | **W:** www.syjlab.com

