

## KASON Series Analytical Balances



### Product instruction

The visual appearance of the product is blended into the simple concept, stylish and generous. Specially makes the fine pouring cast- plastic outer covering, so as to enhance the balance of the anti-static anti-jamming capability The panel with five operational button, Sliding glass windscreen provide big space, insightful vision and operate easily. We use high precision sensor to keep the weighing more fast and stable.

### Technical data:

MODEL	KASON-M103	KASON-M203	KASON-M303	KASON-M603
Capacity(g)	110g	210g	310g	600g
Resolution(g)	0.1g/ 0.01g			
Min weighing(g)	0.04			
LINER	±0.003g			
Stable Time	≤2s			
Operation temp	17.5℃~22.5℃			
Repeat ability	0.02g			
Liner	0.03g			
Option	RS232C & RS485,Printer			

# FOCUS IN MATERIAL TEST

**KASONTTEST**®

JINAN KASON TESTING  
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

**P:** +86 159 1008 1986

**E:** admin@jnkason.com | **W:** www.syjlab.com

