

Mechanical Stability Tester for Natural Rubber Latex Concentrate



Product instruction

This tester supplies a method for the determination of the mechanical stability of rubber latex concentrate, also prevulcanized rubber latex concentrate. It is the essential instrument required to improve the quality of rubber latex.

Technical data:

Parameter	Specification
Display Error	$\leq 0.1\%$
Measurement Resolution	1 r/min
Display Mode	Frequency conversion LED digital display
Standard Test Speed	(14000 ± 200) r/min
Maximum Speed	20000 r/min
Minimum Speed	2000 r/min

FOCUS IN MATERIAL TEST

KASONTTEST®

JINAN KASON TESTING
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | **W:** www.syjlab.com

