

KASON DYQ-1T

Metallographic Specimen Grinding Polishing Liquid Dispenser

DIGITAL DROPLET DISPENSER



★ TECHNICAL PARAMETERS

Model	KASON DYQ-1T
Dispensing speed	About 1 ~ 999 drops / min., can be continuously adjusted
Dispensing time	1-999 minutes, can be continuously adjusted.
Installing and fixing way	Desktop
Power supply	AC 220V, single-phase
Dimensions (L*W*H)	200*170*260mm
Net Weigh	5Kg
Package size	550*480*450mm
Gross weigh	18KG



Metallographic products

Product details



Metallographic consumables



★ PRODUCT FEATURES

- Touch screen settings and automatic control.
- Set and precisely control the droplet speed.
- Set and control the drip time.
- Automatic cleaning dropper.
- It can be connected to a polishing machine and can be automatically controlled.
- One-to-many control by grinding & polishing machine (optional).
- Polishing solution magnetic stirring function, independent speed adjustment (optional).