

KASON P-1

Manual Metallographic Specimen Polishing Machine

POLISHING MACHINE



★ TECHNICAL PARAMETERS

Model	KASON P-1
Disc diameter	Ø203mm
Number of grinding and polishing discs	Single disk
Rotating speed	1400 rpm (can be customized to 900 rpm)
Power supply	180W
Dimensions and Net weight	380*490*350mm, 11kg
Package Size and Gross weight	575*465*540mm, 22.15kg
Power supply	220V, 50HZ, Single phase

★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Metallographic sample polishing machine	1	Product Manual	1
Dust cover	1	Product certification	1
Polished velvet cloth Ø203mm	2	Aluminum mouth ring Ø203mm	1
Drain pipe	1		



Product dimensions

Product details



★ PRODUCT FEATURES

- Speed control system optimized for grinding and polishing machines
- Easy-to-operate switches and knobs for easy operation
- Finely ground and surface-treated work disc to ensure smooth sample surface
- The integrally formed high-strength fiberglass material shell is sturdy and durable and will never rust.



Metallographic consumables