

KASON P-2

Metallographic Specimen Polishing Machine

POLISHING MACHINE



★ TECHNICAL PARAMETERS

Model	KASON P-2
Polishing disc diameter	Ø203mm
Structure	Double disk
Rotating speed	1400rpm (900rpm could be customized)
Motor	180w *2
Dimensions and Net weight	900*480*1027mm,26kg
Package size and Gross weight	1006*565*462mm,41kg
Power supply	AC380v,50Hz,3P

★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Polishing machine working platform	1	Standard accessories	1
Outrigger	2	Dust cover	1
Outrigger beam Ø203mm	1	Flannel Ø203mm	2
Product Manual	1	Selling Ø203mm	2
Product certification	1	Buckle Ø203mm	2



Product details

Metallographic products



Metallographic consumables



★ PRODUCT FEATURES

- Speed control system optimized for grinding and polishing machines
- Easy-to-operate switches and knobs for easy operation
- Finely ground and surface-treated work disc to ensure smooth sample surface
- The integrally formed high-strength fiberglass material shell is sturdy and durable and will never rust.