

KASON P-2G

Metallographic Specimen Polishing Machine



★ TECHNICAL PARAMETERS

Model	KASON P-2G
Polishing disc diameter	Ø203mm
Structure	Double disk
Rotating speed	1400rpm (900rpm could be customized)
Motor	180w *2
Dimensions and Net weight	775*475*952mm,58.9kg
Package size and Gross weight	965*550*1148mm,83.1kg
Voltage	220V, 50HZ, 1 phase

★ STANDARD CONFIGURATION

Name	Qty	Name	Qty
Metallographic sample polishing machine	1	Product Manual	1
Dust cover	1	Product certification	1
Polished velvet cloth Ø203mm	2	Water retaining ring	2
Dust cover	2	Cover button	2
Drain pipe	2	Adjustment screw	4
Water inlet pipe	1		



Metallographic products



Product details



Metallographic consumables

★ PRODUCT FEATURES

- Speed control system optimized for grinding and polishing machines
- Easy-to-operate switches and knobs for easy operation
- Finely ground and surface-treated work disc to ensure smooth sample surface
- The integrally formed high-strength fiberglass material shell is sturdy and durable and will never rust.