

KASON-FPS114 Flat Pressure Sampler

(The picture is for reference only)

Product Introduction

This sampler is an essential sampling tool for the Cobb absorption test. It can quickly and accurately cut samples of a specified area.

Technical Parameters:

Sampling Diameter: 114 mm

Product Features:

Simple instrument structure, high accuracy, and easy maintenance.

FOCUS IN MATERIAL TEST

KASONTTEST®

JINAN KASON TESTING
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | W: www.syjlab.com

