

KASON HTPR-150D Digital Plastic Rockwell Hardness Tester



Large testing table



V-shaped testing table



Medium testing table



Ball indenter



Diamond Rockwell indenter



Standard hardness block

PRODUCT FEATURES:

- Lever Loading, durable and reliable, the test process automation, no human operator error.
- No friction spindle, high precision test force.3.Precision hydraulic buffers, steady load.
- Dial display the hardness value,HRE,HRL,HRM,HRR,and can choose other Rockwell Scale.

TECHNICAL PARAMETERS:

Model	HTPR-150/ XHR-150
Measuring range	70-100HRE, 50-115HRL, 50-115HRR, 50-115HRM
Initial Test Force	98.07N (10Kg)
Total Test Force	588.4, 980.7, 1471N (60, 100, 150kgf)
Max. height of test piece	210mm
Throat	165mm
Min. scale value	0.5HR
Dwelling time	1~99s
Hardness Display	Dial Gauge
Measuring scale	HRG, HRH, HRE, HRK, HRL, HRM, HRP, HRR, HRS, HRV
Power supply	220V AC or 110V AC, 50 or 60Hz
Dimensions	527*242*702mm
Weight	90kg
Package size	620*490*826mm
Gross weight	102kg

STANDARD CONFIGURATION:

Rockwell hardness tester host	1 Set	Power cable	1 Pc
Diamond Rockwell indenter	1 Pc	Dust-proof cover	1 Pc
Ball indenter	3 Pc	Accessory case	1 Pc
Standard hardness block	4 Pcs	Operating instruction	1 Pc
Large testing table	1 Pc	Quality certificate	1 Pc
Medium testing table	1 Pc	Packing List	1 Pc
V-shaped testing table	1 Pc	Standard Weight	3 Pcs