

**KASON SC12 Sample Cutters for Plastic Film****KASONTTEST®****Application**

The Sample Cutter is used for cutting samples of various plastic films, paper and other materials prior to stretching, peeling and heat-sealing tests. It is easy and efficient to operate and can accurately cut samples of standard required sizes.

**Features**

- 1.Meets the sample making requirements for Type 2 specimens of standard GB/T 1040.3.
- 2.Using imported tungsten steel finely ground special blades, sharp and durable, with hardness up to 92 HRA.
- 3.Accurate sample size, no damage, smooth and burr-free edges.
- 4.Spring tensioned for easy sample fixing and smooth sample cutting.

- 5.The ability to cut up to 12 samples in a single pass allows for efficient production of large batches of samples.
- 6.Safety protection, the outer layer of the blade is covered with an aricotta stopper so that the blade does not come into direct contact with it.
- 7.Different sample widths can be ordered to suit the requirements of each customer.

**Technical Parameters**

item	Technical parameters
Cut size	15 mm×210 mm (Optional width)
Number of samples	12 pcs (Single Layer)
Sample thickness	≤ 200 μm
Product size	250 mm ×350 mm ×248 mm

# FOCUS IN MATERIAL TEST

## **KASONTTEST**®

JINAN KASON TESTING  
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | W: www.syjlab.com

