

Universal Wear Abrasion Tester

KASONTTEST®



Note: The picture is for reference only

Application

The universal wear tester is a multi-functional and precise abrasion tester that can adapt to the abrasion resistance of shoelaces, fabrics, PVC, belts and other materials. It can test wear resistance, bending and wear of the edges of the material on the surface of the material. The test stops automatically.

Compliance with the standard: in accordance with ASTM-D3885 AATCC-119&120 ASTM-D3886 ASTM-D3514 ASTM-D3885 DRAFT 20-0-70 standard.

Technical Parameter:

Pressure range	1~9psi
Travel	25 mm \pm 2 mm
Pull the weight	0.5 lb, 1 lb, 2 lb.
Balance weight:	0.5 lb, 1 lb, 2 lb.
Counter	0~99999
Timer	30min
Reciprocating speed:	115 times / min \pm 10
Volume (W*D*H)	31*65*67cm

FOCUS IN MATERIAL TEST

KASONTTEST®

JINAN KASON TESTING
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | **W:** www.syjlab.com

