

KASON-TG4 Electric Thickness Gauge



(The picture is for reference only)

Product Introduction

This electric thickness gauge is suitable for measuring the thickness of flat sheet samples. Because it eliminates human error, it is more suitable for testing the thickness of thin paper, offering better accuracy and repeatability.

Standards: ISO534, GB451.3, QB/T1055 Applicable to paper and paperboard thickness measurement

Technical Parameters:

1. Measuring Range: 0-4 mm
2. Dial Graduation: 0.001 mm
3. Accuracy: ± 0.0025 mm or $\pm 0.5\%$
4. Contact Pressure: 100 ± 10 kPa
5. Contact Area: 2 ± 0.05 cm²

6. Dimensions: 150 x 120 x 350 (mm)
7. Power Supply: 220V
8. Testing Frequency: 9 times/minute
9. Dimensions: 130 × 160 × 360 mm (L x W x H)
10. Net Weight: 15 kg

Product Features:

Digital dial indicator, motorized, capable of continuous automatic testing.

FOCUS IN MATERIAL TEST

KASONTEST®

JINAN KASON TESTING
EQUIPMENT Co, LTD.

DuandianIndustrial Park , Jingshi Road, Jinan City,China.

P: +86 159 1008 1986

E: admin@jnkason.com | **W:** www.syjlab.com

